## A New Standard of Oral Care for the NICU

Plak-Vac Petite Oral Hygiene Kits are tailor-made to support the oral hygiene needs of both intubated and non-intubated babies. Our Petite kits are configured to include choices for swabbing, cleaning, moisturizing, and suctioning the oral cavity. Organizing all components at the point-of-care enhances protocol compliance as well as adherence to infection control policies.









449 Sovereign Court St. Louis, MO 63011

Web: www.trademarkmedical.com Phone: 636-527-2288

Made in U.S.A. M1017H

## **Distributed By:**

CONNECT MEDICAL CORPORATION 118 Homestead Drive Suite 3 Phone: 1-866-416-NICU (6428) • Fax: 1-866-415-8962 www.connectmedical.biz • info@connectmedical.biz

## **Plak-Vac Petite Swab Options** Low Absorption Head 6 mm Head 8 mm Head Actual Head Size 12 mm Head



## Petite Oral Hygiene Kit Worksheet Treatment $oxedsymbol{\square}$ Daily Organized Kit Kit Component Options Quantity/Kit Frequency Plak-Vac® Petite Swabs Treatment Frequency (q) Petite Oral Swab Low Absorption Specific Time (ie 08:00) Petite Oral Swab 6 mm Petite Oral Swab 8 mm Time Range Petite Oral Swab 12 mm (ie 10:00- 12:00) Cleansers/Moisturizers Time Caption 5 ml Sterile Water □ Saline□ (ie Hours, Clean, Suction) Biotene\* Oral Balance Gel 3g Saline Wipes Kit Start Time White Petroleum Jelly 5g Packet Lansinoh\* Lanolin 1.5 g Packet Multi-Day Kit Ayr\* Saline Nasal Gel 0.5 oz. Tube Oral Suction & Other Small Beginnings Bebeonker\* Neotech Little Sucker\* (Specify Other Labeling/Instructions: Suction Catheter 8 / 10 / 12 Fr. 2" by 2" Gauze Medicine Cup Please Submit Protocol to your Sales Representative

Disclaimer: Trademark Medical makes no claim of safety or efficacy for any products used in Plak-Vac Oral Care Kits beyond those made by the product manufacturer's legally approved claims and indications for use.